


Revolutionary.

Introducing the OnControl Bone Marrow System, the first REAL advance for bone marrow procedures in 50 years.

“The OnControl System has provided the Department of Pathology with very generous core biopsies of bone. The samples are the best that Dr. Darling and I have ever seen. In sum, the Department of Pathology is extremely pleased with the quality of specimen that OnControl provides, and we look forward to (rather than dreading) bone biopsies in the future. This device, I guarantee, represents an important improvement in good patient care.”


– Dr. Martinez, Pathologist



The sections shown are from samples of bone marrow biopsy obtained using the OnControl Bone Marrow Biopsy 11g needle. The sample measured 1.7 cm.

Advancing the Process,
Enhancing the OutcomeSM

OnControl™ by **vidacare**[®]
BONE MARROW SYSTEM



OnControl™

For inquiries or more information:

Customer Service1-866-403-1287
Corporate Office1-210-375-8500
EmailVAsales@Vidacare.com
Web Vidacare.com

vidacare[®]


4350 Lockhill Selma Road, Suite 150
Shavano Park, Texas 78249



Use your mobile device to go to our website.
Get the free mobile app at <http://gettag.mobi>

M-401 Rev A

98% successful core capture*



Old Needle Core Sample

OnControl Core Sample

Bone Marrow Biopsy System


Advancing the Process,
Enhancing the OutcomeSM

Vidacare.com

Quality.

Core Capture has a new standard.

After a lithium-powered driver advances the needle with speed, ease and precision into the medullary space, the inner stylet is removed, and the threaded cannula effectively “grabs” the marrow. The result is exceptional quality core samples** – quickly and consistently.



 **nControl™**. Making a powerful difference.

“As a physician who performs [traditional] bone marrow procedures regularly, I know how physically hard it can be. For nearly 3 years, I’ve suffered from pain in my right (dominant) hand, and doing procedures exacerbated that pain. I have had a great relief with my wrist pain using the OnControl system. After months of unhelpful pain management measures, I was able to get my hands on the OnControl system, and it helped immediately, just by taking the strain off of my injured/damaged tendon. Now, just a few months later, the pain is virtually gone, and it has made a huge difference in my life.”

– Dr. McCaul, Oncologist

When you perform a Bone Marrow Biopsy or Aspiration, it’s all about the sample – for you, and most importantly, for your patient.

- Easy-to-Use
- Powerful
- Accurate

* Brenner A. 2010, “Comparing Powered Bone Marrow Biopsy Procedures to Manual Bone Marrow Biopsy Procedures Using Healthy Volunteers,” 52nd ASH Annual Meeting and Exposition, Orlando, FL, Orange County Convention Center, December 4-7, 2010. Research sponsored by Vidacare Corporation.

** Swords A. 2010, “A New Rotary Powered Device for Bone Marrow Aspiration and Biopsy Yields Superior Specimens with Less Pain: Results of a Randomized Clinical Study,” 52nd ASH Annual Meeting and Exposition, Orlando, FL, Orange County Convention Center, December 4-7, 2010. Research sponsored by Vidacare Corporation.

THE ONCONTROL BONE MARROW BIOPSY SYSTEM

