

EZ-IO® Product Systems for Military Use

FEATURES

- Seven FDA cleared sites (sternum, humerus, tibia)
- Battery powered or manually operated
- Powered insertion under 10 seconds
- Manual insertion under 30 seconds
- Standard Luer Lock
- Drug delivery equivalent to central venous access
- Simple removal – no special tools needed
- Complete realistic training (mannequin, CD's, and website)

SPECIFICATIONS

EZ-IO G3 Power Driver – Sealed Lithium Battery (Reusable) (Not for sternal use)

Power	Sealed lithium (no battery change required)
Dimensions	16.5 cm x 11.4 cm x 6.4 cm (6.5 in x 4.5 in x 2.5 in)
Weight	315 grams (11.1 ounces)

EZ-IO 15 mm Needle Set (Disposable) (Not for sternal use)

Size	15 gauge (1.8 mm) x 15 mm long
Weight	5 grams (in sterile pack)
Material	304 stainless steel
Other	Latex free

EZ-IO 25 mm Needle Set (Disposable) (Not for sternal use)

Size	15 gauge (1.8 mm) x 25 mm long
Weight	5 grams (in sterile pack)
Material	304 stainless steel
Other	Latex free

EZ-IO 45 mm Needle Set (Disposable) (Not for sternal use)

Size	15 gauge (1.8 mm) x 45 mm long
Weight	5 grams
Material	304 stainless steel
Other	Latex free

Manual Needle Set Kit (Not for sternal use)

Kit Dimensions	15.2 cm x 8.3 cm (6 in x 3.25 in) foldable
Kit Weight	40 grams (1.4 ounces)
Size	15 gauge (1.8 mm) x 25 mm long
Material	304 stainless steel
Other	Latex free

Sternal Needle Set Kit

Kit Dimensions	15.2 cm x 8.3 cm (6 in x 3.25 in) foldable
Kit Weight	45 grams (1.6 ounces)
Size	15 gauge (1.8 mm) x 7.5 mm long
Material	304 stainless steel
Other	Latex free




PROXIMAL TIBIA



PROXIMAL HUMERUS



DISTAL TIBIA



STERNUM

EZ-IO® Product System Ordering Information

Part #	NSN	Description	Contents
9050	6515-01-537-9615	EZ-IO Power Driver (no case)	1 Driver
9040	6515-01-571-3152	G3 Power Driver (no case)	1 Driver
9001	6515-01-537-9007	25 mm Needle Sets - for patients 40 kg and greater. Size: 25 mm; 15 gauge	5 25 mm Sterile Needle Sets 5 EZ-Connects 5 Wrist Bands 1 Directions for Use
9018	6515-01-537-9013	15 mm Needle Sets - for patients 3 to 39 kg. Size: 15 mm; 15 gauge	5 15 mm Sterile Needle Sets 5 EZ-Connects 5 Wrist Bands 1 Directions for Use
9026	6515-01-541-1946	Military Sterile Training Needle Sets	5 25 mm Sterile Training Needles/Box
9079	6515-01-577-0312	45 mm Needle Sets - for patients 40 kg and greater. Size: 45 mm; 15 gauge	5 EZ-Connects 5 45 mm Sterile Needle Sets 1 Directions for Use
9017	6515-01-540-9808	25 mm Training Needle Sets - for patients 40 kg and greater	1 25 mm Training Needle Set
9021	6515-01-541-9382	15 mm Training Needle Sets - for patients 3 to 39 kg	1 15 mm Training Needle Set
9078	Pending NSN	45 mm Training Needle - for patients 40 kg and greater	1 45 mm Training Needle Set
9045	6515-01-559-7489	Sternal Needle Set Kits	1 Sternal Needle Set 2 Alcohol Preps 1 2x2 1 Wrist Band 1 Lancet 1 Directions for Use 1 Sharps Protector
9047	6910-01-559-7065	Sternal Training Needle Set Kits	1 Sternal Needle Set 2 Alcohol Preps 1 EZ-Connect 1 2x2 1 Wrist Band 1 Lancet 1 Directions for Use 1 Sharps Protector
9055	6515-01-557-6937	Manual Needle Set Kits	1 Manual Needle Set 2 Alcohol Preps 1 EZ-Connect 1 2x2 1 Wrist Band 1 Directions for Use 1 Sharps Protector
9057	6515-01-557-6962	Manual Training Needle Set Kits	1 Manual Needle Set 2 Alcohol Preps 1 2x2 1 Sharps Protector 1 EZ-Connect 1 Directions for Use
9034		G3 Training Kit complete power IO Training Kit	1 G3 Training Driver 3 Needle Sets and Bones
9044	Pending NSN	Sternal Training Mannequin	1 Sternal Training Mannequin
9042	Pending NSN	Sternal Training Replacement Bones and Skin	2 Bones, 1 Skin
9037	Pending NSN	Soft Tissue IO Extremity Trainers (Adult Tibia, Pediatric Tibia)	2 Lower Leg Mannequins
9041	Pending NSN	Training Bones (Adult Proximal Tibia, Distal Tibia, Proximal Humerus, Pediatric Tibia)	1 Adult Proximal Tibia 1 Distal Tibia 1 Proximal Humerus 1 Pediatric Tibia

ACCESSORIES

9065	Pending NSN	G3 Vascular Access Pack (9040 only)	1 Vascular Access Pack
------	-------------	-------------------------------------	------------------------

ORDERING THE EZ-IO: Vidacare.com military@vidacare.com

For military orders outside the US, consult our website (Vidacare.com) for your local Vidacare distributor.