

The OnControl™ Bone Access System

FEATURES

- Easy to use
- Improved tactile feedbackⁱ
- Fast insertion
- Battery powered insertion
- Standard Luer-lock connection
- Effective and safe

SPECIFICATIONS

OnControl Power Driver – Sealed Lithium Battery (Reusable)

Power	Sealed Lithium (no battery change required) Battery Indicator Light
Dimensions	16.5cm x 11.4cm x 6.4cm (6.5in x 4.5in x 2.5in)
Weight	315g (11.1oz)

OnControl Bone Access Needle Set 152mm

Size	152mm (6in)
Gauge	11ga
Outer Diameter	3.048mm (0.12in)
Inner Diameter	2.54mm (0.1in)
Material	304 Stainless Steel
Other	Latex Free

Sterile Procedure Tray

Dimensions	31.75cm x 25.4cm x 5.7cm (12.5in x 10in x 2.25in)
Weight	289g (10.2oz)
Other	Latex Free

REFERENCES

- ⁱ Miller L. 2010, "Tactile Feedback Comparison of Manual, Hammer, and Rotary Powered Needle Insertion Accuracy for Bone Marrow Biopsy and Other Procedures," 52nd ASH Annual Meeting and Exposition, Orlando, FL, Orange County Convention Center, December 4-7, 2010. Research sponsored by Vidacare Corporation.

