

EZ-IO® Intraosseous Infusion System

FEATURES

- For patients 3kg and larger
- Insertion under 10 secondsⁱ
- Stable and secure
- Drug delivery equivalent to central linesⁱⁱ
- Effective, safe, multi-site placement
- Battery power makes insertion easy and controlled
- Standard Luer-lock catheter
- Simple removal

SPECIFICATIONS

EZ-IO G3 Power Driver – Sealed Lithium Battery (Reusable) – NEW

Power	Sealed Lithium (no battery change required)
Dimensions	16.5cm x 11.4cm x 6.4cm (6.5in x 4.5in x 2.5in)
Weight	315 grams (11.1 ounces)

EZ-IO Power Driver (9050) – Sealed Lithium Battery (Reusable)

Power	Sealed Lithium (no battery change required)
Dimensions	23cm x 13cm x 8cm (9in x 5in x 3in)
Weight	456 grams (16.1 ounces)

EZ-IO 15mm Needle Set (Disposable)

Size	15 gauge (1.8mm) x 15mm long
Weight	5 grams
Material	304 Stainless Steel
Other	Latex Free

EZ-IO 25mm Needle Set (Disposable)


Size	15 gauge (1.8mm) x 25mm long
Weight	5 grams
Material	304 Stainless Steel
Other	Latex Free

EZ-IO 45mm Needle Set (Disposable) – NEW

Size	15 gauge (1.8mm) x 45mm long
Weight	5 grams
Material	304 Stainless Steel
Other	Latex Free

EZ-Stabilizer

Dimensions	14.6cm x 13.3cm x 3.2cm (5.75in x 5.25in x 1.25in)
Weight	17 grams (0.602 ounces)
Other	Latex Free



REFERENCES

ⁱ Cooper BR, Mahoney PF, Hodgetts TJ, Mellor A. Intra-osseous access (EZ-IO®) for resuscitation: UK military combat experience. *JR Army Med Corps.* 2008;153:314-6.

ⁱⁱ Hoskins S, Kramer G, Stephens C, Zachariah B. Efficacy of epinephrine delivery via the intraosseous humeral head route during CPR. *Circulation.* 2006;114:II_1204.

EZ-IO® Product System Ordering Information

Part #	Description	Contents
9001	25mm Needle Sets – for patients 40kg and greater Size: 25mm; 15 gauge	5 25mm Sterile Needle Sets 5 EZ-Connects 5 Wrist Bands 1 Directions for Use
9018	15mm Needle Sets – for patients 3 to 39kg Size: 15mm; 15 gauge	5 15mm Sterile Needle Sets 5 EZ-Connects 5 Wrist Bands 1 Directions for Use
9079	45mm Needle Sets – for patients 40kg and greater Size: 45mm; 15 gauge	5 45mm Sterile Needle Sets 5 EZ-Connects 5 Wrist Bands 1 Directions for Use
9017	25mm Training Needle	1 25mm Training Needle Set
9021	15mm Training Needle	1 15mm Training Needle Set
9078	45mm Training Needle	1 45mm Training Needle Set
9058	G3 Power Driver – Sealed Lithium Battery	1 Driver
9059	G3 Power Training Driver – Not for Clinical Use	1 Driver
9050	Power Driver – Sealed Lithium Battery	1 Driver
9052	Power Training Driver – Not for Clinical Use	1 Driver

ACCESSORIES

G3 Part #	Description	Contents
9072	G3 Hard-sided Carrying Case (9058 & 9059 only)	1 Case
9074	G3 Crash Cart or wall/bag Mountable Driver Cradle (9058 & 9059 only)	1 Cradle
9065	G3 Vascular Access Pack (9058 & 9059 only)	1 Vascular Access Pack



9050 Part #	Description	Contents
9007	Hard-sided Carrying Case (9050 & 9052 only)	1 Case
9029	Crash Cart or wall/bag Mountable Driver Cradle (9050 & 9052 only)	1 Cradle
9027	Vascular Access Pack (9050 & 9052 only)	1 Vascular Access Pack



Products Part #	Description	Contents
9013	EZ-Connect 90° Low Profile Extension Sets	50 EZ-Connect Sets
9066	EZ-Stabilizer	1 EZ-Stabilizer

